



CERTIFICATE OF ANALYSIS

Product PEROXIDASE WHEAT GERM AGGLUTININ

Catalog No. PL-1026

Amount 2 mg active conjugate

Lot No. W1026

Concentration 40 mg/ml

Solution buffered salt, 0.02% Thimerosal

Active Conjugate (as protein) 98%

Molar ratio of active conjugate $\frac{\text{Moles of HRP}}{\text{Moles of WGA}} = \underline{1.0}$

Storage Conditions Refrigerate

Remarks:

This lectin conjugate has been specially prepared for neuronal transport studies where a high concentration of protein is required.