

CERTIFICATE OF ANALYSIS

Product RHODAMINE CONCANAVALIN A

****AFFINITY PURIFIED****

Catalog No. RL-1002

Amount 25 mg protein

Lot No. <u>ZB1124</u>

Protein Concentration 5.0 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.1 mM Ca++

0.08% sodium azide, 0.01 mM Mn++

Free Rhodamine None

F/P (molar) 3.4

Storage Refrigerate in the dark

Other Remarks:

If precipitate forms upon long term storage at 4 °C, warm to 37 °C and centrifuge before use. This does not have a significant adverse effect on the product.

The recommended concentration range for use is 5-20 µg/ml.