



---

## CERTIFICATE OF ANALYSIS

Product **RHODAMINE DOLICHOS BIFLORUS AGGLUTININ**

Catalog No. RL-1032

Amount 2 mg

Lot No. ZB0821

Protein Concentration 2 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide, 0.1 mM Ca<sup>++</sup>

Free Rhodamine None

F/P (molar) 3.8

Maximum Excitation @ 545-555nm

Maximum Emission @ 570-580 nm

Storage Refrigerate in the dark

Other Remarks:

If precipitate forms upon long term storage, warm to 37 °C.

The recommended concentration range for use is 5-20 µg/ml.