

CERTIFICATE OF ANALYSIS

Product RHODAMINE ULEX EUROPAEUS AGGLUTININ I

Catalog No. RL-1062

Amount 2 mg

Lot No. ZB1001

Protein Concentration 2.0 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,

Free Rhodamine None

F/P (molar) 3.2

Storage Refrigerate in the dark

Other Remarks:

A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage at 4 °C, warm to 37 °C and centrifuge before use.