



CERTIFICATE OF ANALYSIS

Product **RHODAMINE PEANUT AGGLUTININ**

Catalog No. RL-1072

Amount 5 mg

Lot No. ZC0606

Protein Concentration 5 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Rhodamine None

F/P (molar) 4.0

Maximum Excitation @ 545-555nm

Maximum Emission @ 570-580 nm

Storage Conditions Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 µg/ml.

Inhibiting/Eluting Sugar: 200 mM galactose (S-9003).

For immunofluorescence applications, we recommend that sections be mounted in VECTASHIELD[®] Hard♦Set[™] Mounting Medium (H-1400) or VECTASHIELD[®] Hard♦Set[™] Mounting Medium with DAPI (H-1500).