



CERTIFICATE OF ANALYSIS

Product **RHODAMINE GRIFONIA (BANDEIRAEA) SIMPLICIFOLIA
LECTIN I (GSL I, BSL I)**

Catalog No. RL-1102

Amount 2 mg

Lot No. ZC0524

Protein Concentration 2.0 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide, 0.1 mM Ca⁺⁺

Free Rhodamine None

F/P (molar) 10.8

Maximum Excitation @ 545-555 nm

Maximum Emission @ 570-580 nm

Storage Conditions Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 µg/ml.