



CERTIFICATE OF ANALYSIS

Product **RHODAMINE -*PHASEOLUS VULGARIS* LEUCOAGGLUTININ (PHA-L)**

Catalog No. RL-1112

Amount 2 mg protein

Lot No. ZA0317

Protein Concentration 2.0 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide, 0.1 mM Ca⁺⁺

Free Rhodamine None

F/P (molar) 3.6

Storage Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 μ g/ml.