



CERTIFICATE OF ANALYSIS

Product **CHITIN HYDROLYSATE**
Catalog No. SP-0090
Amount 10 ml
Lot No. U1009
Solution Nearly saturated sodium chloride in water
Storage Conditions Refrigerate

Remarks

This lot of Chitin Hydrolysate contains a nondialyzable O.D. _{280nm} of no greater than 0.1.

If a precipitate forms during storage do not warm solution. This may cause further precipitation. Centrifuge or filter the hydrolysate to remove precipitated material.

For additional instructions on use see other side.