



CERTIFICATE OF ANALYSIS

Product **ANTI-AVIDIN D**
** AFFINITY PURIFIED **

Catalog No. SP-2000

Amount 1 mg

Lot No. V0705

Produced in goat

Solution Reconstitute by adding 1 ml of water. The resulting solution will have the following composition: 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage Refrigerate; for long-term storage, aliquot and store frozen.

Antibody activity by solid phase binding assay 97%

Remarks:

The recommended concentration range for use is 2-5 µg/ml.