



CERTIFICATE OF ANALYSIS

Product ANTI-BIOTIN
** AFFINITY PURIFIED **

Catalog No. SP-3000

Amount 1 mg

Lot No. Y2504

Produced in goat

Solution Reconstitute by adding 1 ml water. The resulting solution will have the following composition: 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage Refrigerate

Antibody activity by solid phase binding assay 99%

Remarks:

This product is recommended for use in transfer blot applications (northern, Southern, western, dot blots, etc.), immunohistochemical, and *in situ* hybridization applications.

The recommended concentration range for use is 0.5-5 µg/ml. This reagent can be detected with a VECTASTAIN[®] ABC System or the ImmPRESS[™] Anti-Goat IgG Polymer Detection System (Cat. No. MP-7405).