



CERTIFICATE OF ANALYSIS

Product FLUORESCCEIN ANTI-BIOTIN
** AFFINITY PURIFIED **

Catalog No. SP-3040

Amount 0.5 mg

Lot No. ZB0611

Produced in goat

Concentration 1 mg/ml

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide

Free Fluorescein none

F/P (molar) 8.3

Storage Refrigerate in the dark

Antibody activity by solid phase binding assay 99%

Remarks:

The recommended concentration range for use is 5-20 µg/ml.