



CERTIFICATE OF ANALYSIS

Product ANTI-STREPTAVIDIN
** AFFINITY PURIFIED **

Catalog No. SP-4000

Amount 1 mg

Concentration 2 mg/ml

Lot No. ZB0710

Produced in goat

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage Refrigerate; for long-term storage, aliquots may be stored frozen.

Antibody activity by solid phase binding assay 98%

Remarks:

The recommended concentration range for use is 2-5 μ g/ml.