



CERTIFICATE OF ANALYSIS

Product TEXAS RED[®] LYCOPERSICON ESCULENTUM
(TOMATO) LECTIN

Catalog No. TL-1176

Amount 1 mg active conjugate

Lot No. ZB0722

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide
0.1 mM Ca⁺⁺.

Concentration 1 mg/ml

Free Fluorochrome none

F/P (molar) 3*

Active conjugate 99%

Storage Refrigerate in the dark.

Other remarks:

If a precipitate forms upon long-term storage, warm to 37 °C.

*J Immunological Methods 50, 193 (1982)

® Trademark of Molecular Probes